



## Rental Agreement

**By this agreement** the owner lets and the borrower takes the equipment on hire in accordance with the attached terms and conditions.

Hirer Details	
Name	
Address	
Contact Number	
Proof of address (Utility Bill or Bank Statement no older than 3 months old)	
Photo ID (Passport or Driving Licence)	
Proof of eligibility (PIP/Motability letter or signed letter from a Medical Professional)	
Registered Hiccups family date	

Hire Period (not longer than 21 days)	
Start date	
Finish date	
Agreed Return Time	





Delta Buggy size (please refer to website for size guidance: )	
Medium (Delta Buggy M2 Navy)	
Large (Delta Buggy L1 Navy)	

Additional items hired				
	Checking Out		Checking In	
	Comments	Initials	Comments	Initials
Delta Frame Side Pad Set Navy				
Delta Strap Pad Set Navy <i>(Size M2 or L1/L2)</i>				
Spare Canvas with Pads & Straps - Black/Orange <i>(Size M2 or L1/L2)</i>				
Delta Sun Canopy Black <i>(Size M2 or L1/L2)</i>				
Delta Rain Cover <i>(Size M2 or L1/L2)</i>				
M2/L/XL Headrest Cushion for Delta Buggy				
M2/L1 Seat Infil Cushion for Delta Buggy				
200mm Lateral Rolls for Delta Buggy				
Pommel for Delta Buggy <i>(Size M2 or L/XL)</i>				





Seat length extender for Delta Buggy (Size M2 or L/XL)				
Delta Buggy Travel Bag (Medium or Large)				
Beach Wheels for Delta Buggy				
Delta Trail-M1/M2 : Cycle Trailer Conversion Kit (Size M1/M2 or L1/L2)				

Delta All Terrain Buggy Hire safety checklists				
	Checking Out		Checking In	
	Comments	Initials	Comments	Initials
Anti-tipper function and secure locking into place; when deployed the end of anti-tipper should be approx. 22mm from the ground.				
Security of the quick release wheel attachments				
Wear to folding joints. If there is excessive movement in the joints, take out a stainless steel screw to inspect the polymer bearing that protects the aluminium joint from the fixing				





screw. The bearings can be replaced if necessary.				
Wear to the handle adjusters clutch mechanism. When the handle does not lock securely the adjusters should be replaced.				
That all screws and fixings are secure; folding joint fixings should not be over tightened.				
Condition and secure attachment of the foam that covers the handle to backrest frame joint.				
Tyre wear and inflation pressure. Standard tyres should be inflated to 20 and 30 p.s.i. (1.4 and 2 bar). Beach wheel tyres should be inflated to between 10 and 15 p.s.i. (0.7 and 1 bar approx).				
Secure attachment of footrest; adjust height if required.				
Check condition of harness and foot-straps; replace if necessary.				
Any other scuff marks / damage				





Instruction manual included				
Terms and Conditions are included				
First Aid kit included				
Accident and Fault Reporting Form included				
Lock included				
Slime and Pump included				





Little Hiccups is a small registered charity and any donation you may wish to make from your deposit would be gratefully received. A donations help us to continue to provide support for children with additional needs. If you would like to donate your deposit, or part of it, please specify below. Please leave blank if not applicable.

Amount (£)	
------------	--

I agree that I am happy to borrow out the Delta All Terrain Buggy and that I have read the Terms and Conditions, the Instruction Manual and have reviewed the condition of the Delta with the Little Hiccups Agent. I can confirm that I have also been shown how to use the Delta safely and will conform to the safe system of work.

Signed out on behalf of Little Hiccups \_\_\_\_\_

Signed out by/on behalf of the borrower \_\_\_\_\_

Date: \_\_\_\_\_

Signed in on behalf of Little Hiccups \_\_\_\_\_

Signed in by/on behalf of the borrower \_\_\_\_\_

Date: \_\_\_\_\_

All information and photos will be stored according to the Data Protection Act 1998 and the GDPR Act 2018. Information collected will be used for Little Hiccups purposes only. We will not share your information with any other companies or organisations.

Please contact us if you have any questions about our privacy notice or information we hold about you. You can contact us by email at [info@littlehiccups.co.uk](mailto:info@littlehiccups.co.uk)

You can see the full privacy notice on our website at <http://www.littlehiccups.co.uk/privacy.htm>

